

FULLY AUTOMATIC AUTOCLAVES TABLE TOP SA-230FA / 260FA SERIES



FEATURES :

- Innovative and versatile autoclave, which combines compact and elegant design with high performance, reliability and safety
- Microprocessor control system. Fully Automatic control from add-water to sterilization, to dry and 9 programs for your professional selection
- Special Liquid Sterilization program, it is the special function for laboratory and scientific use.

SPECIFICATIONS:

MODEL :	SA-230FA (Class S)	SA-260FA (Class S)
Overall Dimensions	594mm(D) x 522mm(W) x 410mm(H)	655mm(D) x 533(W) x 442mm(H)
Chamber Dimension	230mm x 410mm	260mm x 450mm
Capacity	16l	24l
Net Weight	48kg	51kg
Program Display	Yes	Yes
Power Supply	230V 50 / 60 Hz	230V 50 / 60 Hz
Power Consumption	1.9KW, 9 Amp	2.7KW, 12 Amp
Auto-fill Water System	Yes	Yes
Water capacity of water tank	4200cc	4200cc
Auto Drying	Yes	Yes
Door Closed, Low Water, Over Heat Display	Yes	Yes
Safety Devices	<ul style="list-style-type: none"> • Pressure Safety Valve • Over Temperature Protection • Over Pressure Protection Switch • Door Lock Indication • Pressure Door Auto Lock • Water Level Sensor (Tank & Chamber) • Electric overload protection (20Amp. No fuse Breaker x 2) 	<ul style="list-style-type: none"> • Pressure Safety Valve • Over Temperature Protection • Over Pressure Protection Switch • Door lock Indication • Pressure Door Auto Lock • Water Level Sensor (Tanks & Chamber) • Electric overload protection (20Amp. No fuse Breaker x 2)
Sterilization Temperature	121°C / 135°C	121°C / 135°C
121°C/135°C, rapped/Unwrapped, dry/undry, Liquid, program	Yes	Yes
Standard Accessories	Sterilization box x 1 ; Heater cover x 1	Sterilization box x 1 ; Heater cover x 1

Optional Accessories	<ul style="list-style-type: none"> • Sterlization Tray Set • Tray x 3 • Tray frame x 1 • Tray holder x 1 • Spring Holder x 2 	<ul style="list-style-type: none"> • Sterlization Tray Set • Tray x 3 • Tray frame x 1 • Tray holder x 1 • Spring Holder x 2
Quality Standard	CE	