

FULLY AUTOMATIC AUTOCLAVES TABLE TOP SA-260MA SERIES



CODE : ASA260MA

FEATURES :

- Microprocessor control system: Fully automatic control from Pre-vacuum to dry and 24 programs for your selection
- Pre-heating by dry heater ensuring chamber temperature balance and perfect drying function
- Special LIQUID sterilization program: it is the special function for laboratory and scientific use
- Special PRION sterilization program: it is the special program for PRION sterilization
- Special BOWIE-DICK test program: the users can identify the perfect sterilization function with this program
- An "ERROR" indicator, to help the users and engineers to find the function problem automatically

SPECIFICATIONS:

MODEL :	SA-260MA
Overall Dimensions	655mm x 533mm x 442mm
Chamber Dimension	260mm x 450mm
Capacity	24 l
Temperature / Pressure	Unwrapped / Wrapped : 121 – 134°C ; 2.2kg per cm ² Liquid : 121°C ; 1.2kg per cm ²
Control System	Pressure control
Power Consumption	220 V: 1.8kW ; 8.3Amp / 240 V: 2.1kW ; 8.8Amp
Auto-fill Water System	Yes
Double Water Level Check Device	Yes : One for water tank ; one for chamber water level check
Auto Drying	Yes : 0, 15, 30, 35 or 40 minutes
Pre-Vacuum	Yes
Dry-Vacuum	Yes
Indicators / Temperature Display	Yes / LED
Sterilization Temp / Time	Wrapped : 121°C - 30 min; 134°C - 15 min Unwrapped : 121°C - 15 min; 134°C - 4 min
Standard Accessories	Sterilization box x 1 ; Tray x 1
Optional Accessories	Tray frame x 1 ; Tray holder x 1 ; Tray x 3 ; Spring x 2
Safety Devices	<ul style="list-style-type: none"> ✓ Temperature control switch ✓ Pressure door auto-lock device ✓ Emergency key – cuts power supply & exhausts chamber pressure ✓ Chamber pressure gauge + pressure & steam pressure safety valve ✓ Door close indication ✓ Low water alarm and indication ✓ Over-heat protector ✓ Safety valve ✓ Electronic circuit safety system ✓ Low water indicator & shutout device ✓ Sterilization completion indicator
Quality Standard	CE-MDD