

ANEROID 5 CUFF SET WITH SINGLE HANDED WITH CHROME GAUGE



CODE : *BPMAP5CS*

Specifications:

- Palm type blood pressure monitor
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Chrome casing 65mm diameter dial
- Black powder-coated handle
- 5 Navy blue washable cuffs. Sizes: Infant, Child, Adult, Obese and Thigh
- Single Tube latex bladder in each cuff
- Standard latex bulb and valve

Packing Details:

Single Unit – 1 complete set in cardboard box
(301mm x 220mm x 700mm), 1.03kg

Carton – Packed per order

Quality Approvals:

CE
ISO