

ANEROID BLOOD PRESSURE METER WITH SINGLE HEAD STETHOSCOPE



CODE : BPMA/STH

Specifications:

- Aneroid and stethoscope combined
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Red zone on gauge between 80mmHg and 160 mmHg
- Navy blue washable adult size nylon cuff with single head stethoscope attached
- 2-tube adult size cuff and bladder
- Standard latex bulb and valve
- Vinyl zipper case

Packing Details:

Single Unit – 1 complete set in cardboard box
(195mm x 110mm x 75mm), 550 grams

Carton – 50 Pieces (550mm x 400mm x 400mm), 25kg

Quality Approvals:

CE
ISO