

DELUXE ANEROID POCKET - GAUGE ON CUFF



CODE : *BPMAPKD*

Specifications:

- Economy type blood pressure monitor
- Shock-proof chrome plated gauge with protective ring
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Red zone on gauge between 80mmHg and 160 mmHg
- Navy blue washable adult size nylon cuff with leather attachment
- 2-tube adult size cuff and bladder
- Standard latex bulb and inflation valve with spring
- Vinyl zipper case

Packing Details:

Single Unit – 1 complete set in cardboard box

(195mm x 110mm x 75mm), 360 grams

Carton – 50 Pieces (560mm x 340mm x 345mm), 23kg

Quality Approvals:

CE

ISO