

# ANEROID BLOOD PRESSURE METER ON MOBILE STAND



**CODE :** *BPMAM*

## **Specifications:**

- Square shaped 14cm diameter dial
- 360° Manometer swivel for easy-reading
- Left and right adjustable
- Aluminum stand and tube
- 5 wheel base with brake
- Adjustable height (range of 850mm to 1250mm)
- Plastic basket at rear of unit
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Red zone on gauge between 80mmHg and 160 mmHg
- Navy blue washable adult size nylon cuff with D-ring
- 2-tube adult size cuff and bladder
- PVC coil tubing
- Standard latex bulb and inflation valve with spring

## **Packing Details:**

**Single Unit** – 1 complete set in cardboard box  
(570mm x 490mm x 160mm), 5kg

**Carton** – 1 Piece (570mm x 490mm x 160mm), 5kg

## **Quality Approvals:**

CE

ISO