


ANEROID PALM SATIN ELITE BLOOD PRESSURE METER



CODE : *BPMAPESCR*

Specifications:

- Satin chrome case and handle, 65mm diameter dial
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Red zone on gauge between 80mmHg and 160mmHg
- Navy blue washable adult size nylon cuff
- Single-tube adult size cuff and bladder
- Unique push button air release valve offers precise mini-deflation
- PVC ribbed bulb – Non Latex
- Large size metal end valve with dust filter
- Vinyl zipper case

Packing Details:

Single Unit – 1 complete set in cardboard box

(190mm x 100mm x 75mm), 550 grams

Carton – 30 Pieces (440mm x 480mm x 325mm), 21.5kg

Quality Approvals:

CE, ISO