

# ANEROID SINGLE HANDED BLOOD PRESSURE METER WITH ABS GAUGE (RED OR BLUE)



**CODE :** *BPMAPABSBL*

### **Specifications:**

- Palm type blood pressure monitor- single handed
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Blue and Red powder-coated handles
- Navy blue or red washable adult size nylon cuff
- 2-tube adult size cuff and bladder
- Standard latex bulb and valve
- Vinyl zipper case- colour co-ordinated

### **Packing Details:**

**Single Unit** – 1 complete set in cardboard box  
(195mm x 110mm x 75mm), 400grams  
**Carton** – 50 Pieces (580mm x 390mm x 400mm), 25kg

### **Quality Approvals:**

CE  
ISO