

# ANEROID SINGLE HANDED PALM DELUXE BLOOD PRESSURE METER WITH CHROME PLATED GAUGE



**CODE :** *BPMASH*

**Specifications:**

- Palm type blood pressure monitor
- Chrome plated metal gauge with 50mm diameter dial  
Aneroid gauge calibrated from 0mmHg to 300mmHg
- Stainless steel handle
- Navy blue washable adult size nylon cuff
- 1 or 2-tube adult size cuff and bladder
- Standard latex bulb and valve
- Vinyl zipper case

**Packing Details:**

- Single Unit** – 1 complete set in cardboard box  
(195mm x 110mm x 75mm), 380 grams
- Carton** – 50 Pieces (580mm x 390mm x 400mm), 23kg

**Quality Approvals:**

CE  
ISO