

ANEROID SQUARE FACE WALL MOUNTABLE BLOOD PRESSURE METER



CODE : *BPMAWMSQ*

Specifications:

- Round shaped 140mm diameter dial
- 90° Manometer swivel for easy-reading
- Left and right adjustable
- Metal basket behind unit
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Navy blue washable adult size nylon cuff with D-ring
- 2-tube adult size cuff and bladder
- PVC coil tubing
- Standard latex bulb and inflation valve with spring

Packing Details:

Single Unit – 1 complete set in cardboard box

(175mm x 175mm x 135mm), 850g

Carton – 20 Pieces (370mm x 890mm x 290mm), 17.5kg

Quality Approvals:

CE
ISO