


BP METER ANEROID WALL SQUARE FACE SWIVEL



CODE : BPAWMSQSW

Specifications:

- Square shaped 14cm diameter dial
- Manometer with horizontal swivel
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Navy blue washable adult size nylon cuff with D-ring
- 2-tube adult size cuff and bladder
- Standard latex bulb and inflation valve with spring
- PVC coiled tubing

Packing Details:

Single Unit – 1 complete set in cardboard box
(165mm x 195mm x 175mm), 960g

Quality Approvals:

CE
ISO