

DELUXE ANEROID WALL/DESK BLOOD PRESSURE METER WITH SQUARE FACE



CODE : *BPMAWDC*

Specifications:

- Square shaped 145 mm diameter dial
- Front and rear adjustable for easy reading
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Navy blue washable adult size nylon cuff with D-ring
- 2-tube adult size cuff and bladder
- Standard latex bulb and inflation valve with spring
- PVC coiled tubing

Packing Details:

Single Unit – 1 complete set in cardboard box

(165mm x 195mm x 175mm), 960g

Carton – 10 Pieces (865mm x 390mm x 405mm), 12kg

Quality Approvals:

CE

ISO