

# STANDARD ANEROID POCKET BLOOD PRESSURE METER



**CODE : *BPMAPKS***

**Specifications:**

- Economy type blood pressure monitor
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Red zone on gauge between 80mmHg and 160mmHg
- Navy blue washable adult size velcro nylon cuff
- 2-tube adult size cuff and bladder
- Standard latex bulb and valve
- Vinyl zipper case

**Packing Details:**

**Single Unit** – 1 complete set in cardboard box

(195mm x 110mm x 75mm), 550 grams

**Carton** – 50 Pieces (560mm x 340mm x 345mm), 21kg

**Quality Approvals:**

CE

ISO