

BP METER ANEROID DUAL CUFF - CHILD (GIRL)



CODE :

Specifications:

- Economy type blood pressure monitor
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Red zone on gauge between 80mmHg and 160 mmHg
- 2-tube infant and child cuff and bladder
- Standard latex bulb and inflation valve with spring
- Nylon Carry Pouch

Packing Details:

Single Unit – 235mm (L) x 140mm (W) x 50mm (D)

Quality Approvals:

CE
ISO