

## "XL" Deluxe Automatic Blood Pressure Monitor

Model : MW701 



- Real Fuzzy Technology
- Hypertension Risk Indicator
- Irregular Heartbeat Detector (IHB)
- 120 memories for 2-users with date and time indication
- Average of last 3 readings
- Jumbo size "XL" digit
- One-touch automatic operation
- Carrying bag
- Latex-free patented cone cuff
- PC Link (Optional)
- Adaptor socket



BHS A/A graded approved



ESH approved



Real Fuzzy Technology



Hypertension Risk Indicator



Irregular Heartbeat Detector (IHB)



120 memories for 2-users



Average of last 3 readings



**Clinically Validated**

This blood pressure monitor has passed the criteria of the European Society of Hypertension by achieving all requirements and achieved the highest possible grade (A/A) according to the criteria of the British Hypertension Society. This monitor further complies with the standard required by the Association for the Advancement of Medical Instrumentation (AAMI).

### Real uzzy Technology

Real Fuzzy Technology determines ideal cuff pressure based on one's systolic blood pressure and arm size. Users no longer need to pre-set the inflation level before measurement. The technology eliminates arm discomfort through unnecessary high inflation settings and also prevents inaccurate readings caused by erroneous cuff inflation levels.



Cone cuff fits arm circumference:  
24-36cm (9.4" ~14.2")  
Other size cuffs available upon request

- Storage and Transportation Condition  
-10°C ~ 60°C (14°F ~ 140°F)  
10% ~ 90% RH