

DELUXE MERCURIAL DESK BLOOD PRESSURE METER



CODE : *BPMMDD*

Specifications:

- Glass tube (3mm) with graduation scale
- Mercurial gauge calibrated from 0mmHg to 300mmHg
- Heavy chrome-plated stainless steel case
- Navy blue washable adult size nylon cuff with D-ring
- 2-tube adult size cuff and bladder with D-ring
- Standard latex bulb and inflation valve with spring

Packing Details:

Single Unit – 1 complete set in cardboard box
(165mm x 195mm x 175mm), 1.1kg

Carton – 20 Pieces (415mm x 315mm x 360mm), 24kg

Quality Approvals:

CE
ISO