

# DELUXE MERCURIAL WALL BLOOD PRESSURE METER WITH FRONT BASKET



**CODE : *BPMMWM***

**Specifications:**

- Glass tube (3mm) with graduation scale
- Mercurial gauge calibrated from 0mmHg to 300mmHg
- Red zone on gauge between 80mmHg and 160 mmHg
- Navy blue washable adult size nylon cuff with D-ring
- 2-tube adult size cuff and bladder with D-ring
- Standard latex bulb and inflation valve with spring
- PVC coil tubing
- Basket for cuff and bulb

**Packing Details:**

**Single Unit** – 1 complete set in cardboard box  
(435mm x 110mm x 180mm), 1.1kg

**Carton** – 6 Pieces (420mm x 400mm x 480mm), 9kg

**Quality Approvals:**

CE  
ISO