

# MOBILE MERCURIAL BLOOD PRESSURE METER ON MOBILE STAND



**CODE : *BPMMM***

## **Specifications:**

- Glass tube (3mm) with graduation scale
- Mercurial gauge calibrated from 0mmHg to 300mmHg
- Red zone on gauge between 80mmHg and 160 mmHg
- Navy blue washable adult size nylon cuff with D-ring
- 2-tube adult size cuff and bladder
- Standard latex bulb and inflation valve with spring
- PVC coil tubing
- Aluminum stand and tube
- 5 wheel base with brake
- Adjustable height (range of 850mm to 1250mm)

## **Packing Details:**

**Single Unit** – 1 complete set in cardboard box  
(485mm x 210mm x 555mm), 9kg

**Carton** – 1 Pieces (485mm x 210mm x 555mm), 9kg

## **Quality Approvals:**

CE  
ISO