



# Resting 12-Lead ECG



## Resting ECG "CardioScan Resting 12-Lead"

### Software

New Computerized Resting ECG 12-Lead ECG with Interval/Amplitudes/Axis measurements/with or without automatic interpretation developed and made by DM Software, implementing all the advanced holter ECG technology from DM Software, and many unique analysis functions:




<b>Version 4.0 Standard</b>	<ul style="list-style-type: none"> <li>■ 12 Lead Resting ECG with or without auto interpretation</li> <li>■ Vector ECG</li> </ul>
<b>Version 4.0 Premier</b>	<ul style="list-style-type: none"> <li>■ 12 Lead Resting ECG with or without auto interpretation</li> <li>■ Vector ECG</li> <li>■ Timed Vector ECG</li> <li>■ SAECG (Signal Averaged ECG) Ventricular Late Potentials analysis</li> <li>■ QTd analysis</li> <li>■ Hour ECG</li> <li>■ High Frequency ECG (HFECG)</li> <li>■ Frequency Cardiogram (FCG)</li> <li>■ Simple Exercise ECG</li> </ul>


- Sampling screen Colour Customization
- Multi ECG Report Formats
- Standard or thermal paper print, B&W or colour print
- Auto Interpretation On or Off
- Measurement Table
- Tests Comparison
- 90s Resting ECG full disclosure
- Multi Hour ECG Monitoring with holter analysis
- Simple Exercise Stress Test with report
- E-mail , PDF, ASCII, HL7, Dicom, GDT Interfaces
- Windows 2000, XP, Vista Compatible, and Mini Laptops running under windows environment...



ECG Sampling Screen



Averaged ECG Complex



ECG Report 3x4+1

## Hardware



12-Lead Blue-tooth Transmitter: Front & Back



Blue-Tooth Receiver



Patient Cable

### Specifications and Features table

Model	12-Lead ECG Blue Tooth Module	12-Lead ECG Blue Tooth Transmitter
Leads	12-Lead, 10 electrodes	12-Lead, 10 electrodes
Sampling Rate	12 Bits, 4096Hz	12 Bits, 4096Hz
Dimension	Trans. 13.5 x 7.5 x 2.4 cm Receiv. 13.9 x 12 x 3 cm	13.5 x 7.5 x 2.4 cm
Weight	Emet. 110 g récept. 240 g	110 g
SAECG * & VCG *	Yes	Yes
FCG *	Yes	Yes
Power	Trans. 4 x AA batt. (6.0V) Receiv. USB Port (5.0V)	4 x AA batteries (6.0V)
ECG Download	Blue Tooth Transmission 2,4 GHz 100 M	Blue Tooth Transmission 2,4 GHz 100 M
Compatibility	CardioScan Resting 12-Lead CardioScan Stress Test	CardioScan Resting 12-Lead CardioScan Stress Test

## Innovative high performance Sampling devices for Desktop, Laptop, Pocket PC

■ New Generation long distance Blue Tooth module wireless tests in private office, hospital or clinique, with Desktop, Laptop, or Tablet PC!

■ Mini ultra light USB module for use for in office tests, and field tests with USB cable to Laptop, Mini Laptop or Windows mobile Pocket PC

■ New Standalone Blue Tooth Transmitter light and transportable, perfect for field tests with Mini Laptops, Windows Mobile Pocket PC

... for all your consultations needs: **Office, Clinique, Hospital, Visits, Emergency!**