


# ELECTRICAL PLASTER SAW



**CODE: PLA001**

## **FUNCTIONS:**

- Whether you work in a busy plaster room, operating theatre, hospital ward, small doctor's clinic or veterinary practice handling large or small animals, this plaster saw will meet your requirements.
- The powerful 150 Watt motor maintains rapid and efficient cutting performance on either plaster or synthetic cast material.

## **FEATURES:**

- High Power
- Heavy duty built-in overload system
- Maximum durability
- Maintenance free
- Robust construction
- Weight 1.6kg
- Voltage 220V
- TV and radio suppressed
- SABS Approved