

Bedside Lockers



D-6

Beside cabinet stainless steel top & base Epoxy coated body
 castors 75 mm
 towel rails x 2
 BOX PACKED SIZE 440 X 530 X 800 MM
 490mm(L)*400mm(W)*755mm(H)
 Epoxy coated steel main body with stainless steel top board and stainless steel bottom base.
 2. With one drawer on drawer glides
 one door with inner shelf board.



D-1 / D-2

ABS Bedside locker
 on castors
 BOX SIZE 510 X 510 X 800 MM
 475*470*755mm Packing size:
 505*500*785mm
 Main features:
 1. ABS engineering plastic material
 2. With one food board, one drawer, one small door (inner with one shelf board)
 3. Tower shelf in lateral side, fold away



D-10

ABS / METAL Bedside locker
 METAL BODY ON WHEELS
 ABS FRONT
 stainless steel top
 BOX SIZE 510 X 510 X 800 MM
 475*470*755mm Packing size:
 505*500*785mm
 Main features:
 1. ABS engineering plastic material
 2. With one food board, one drawer, one small door (inner with one shelf board)
 3. Tower shelf in lateral side, fold away