

CODE

PRODUCT DESCRIPTION

Hospital bed Electric - ICU / High Care



DA-2

5 Function Electric Bed Complete with fold up sides

Estimated dimensions 2180 x 970 x 540 mm

BED HEIGHT 540 - 770 mm.

High quality cold-rolled steel.

Head and legs power adjustment:

Head :0-75 °, legs :0-40 °

Bed head and bed foot tilted :0-20 °

Folding safety rails.

Central Locking swiveling Castor

Bed head and bed foot ends removable .

Emergency CPR function.

Bed: 1980 * 900mm,

Bed body: 2090 * 910 mm.

Weight Capacity: 400 lbs(182 kg).

Net Weight: 128KG.

Battery backup



Five-function Electric Bed
DA-2





DA-10

Multifunction Function Electric Bed

Complete with fold up sides

Weight capacity 150 kg

Estimated dimensions

2180 x 970 x 570 mm

BED HEIGHT 540 - 770 mm.

High quality cold-rolled steel.

All functions have power adjustments:

Head :0-75 °, legs foot down :0-85 °

Trendelenburg head and foot down

tilting :0-15 °

Knee up angle:0-12 °

Folding safety rails.

Central Locking swiveling Castor

Bed head and bed foot ends removable

Emergency CPR function.

Bed: 1980 * 900mm,

Bed body: 1880 * 910 mm.

Weight Capacity: 400 lbs(182 kg).

Net Weight: 128KG.

Battery backup

Voltage:220V

Current:8 A

