

Hospital bed Electric

Hospital bed Electric - Home care, Frail care, Clinic and General Wards



FS3231WM

Bed with battery backup, sides and over bed pole.

- Wooden headboard and footboard
- Folding sides
- Electric backrest and knee section
- Dual motor elevating mechanism
- FS577 Overbed lifting pole and handle.
- Steel epoxy coated bed frame. with rubber brake castors.
- Options: mattresses.
- Package Size: 200cmX96cmX26cm



FS3231WM

- ELECTRIC FIVE FUNCTION METAL NET PANEL
- WOODEN HEAD AND FOOT BED-BOARDS
- ALUMINUM REMOVABLE GUARDRAILS
- POWDER COATING STEEL FRAME
- PEDAL LOCKING CASTORS
- DOUBLE MOTOR LIFTING SYSTEM
- SIZE: 213X101X34





FS577

FS3237WMF3

Bed with battery backup, sides and over bed pole.

- Folded steel slat mattress platform.
- Wooden headboard and footboards.
- Steel folding bed sides.
- Steel epoxy coated bed frame with rubber brake castors.
- Electric backrest and knee section
- FS577 Overbed lifting pole and handle.
- Dual motor elevating mechanism

- Options: mattresses.
- Package Size: 208cmX96cmX26cm



FS3237WMF3

- ELECTRIC FIVE FUNCTIONS MILD STEEL PANEL
- WOODEN HEAD & FOOT BED-BOARDS
- MILD STEEL REMOVABLE GUARDRAILS
- POWDER COATING MILD STEEL FRAME
- PEDAL LOCKING CASTORS
- SIZE :209X100X28





FS3230W

Bed with battery backup, sides and over bed pole.

- Aluminum folding sides.
- Steel epoxy coated bed frame with rubber brake castors.
- Electric backrest and knee section.
- Electric height adjustment.
- FS577 Overbed lifting pole and handle.

- Options: mattresses,
- Package Size: 205cmX96cmX45cm: 103cmX54cmX13cm



FS577

> FS3230W

- ELECTRIC THREE FUNCTIONS, METAL NET PANEL
- ALUMINUM REMOVABLE GUARDRAILS
- POWDER COATING STEEL FRAME
- PEDAL LOCKING CASTORS
- ELECTRIC LIFTING SYSTEM WITH CE,UL APPROVAL
- SIZE: 209X102X49





FS3230WMGF3

Bed with battery backup, sides and over bed pole.

- Folded steel mattress platform.
- Wooden headboard and footboards.
- Steel folding bed sides.
- Steel epoxy coated bed frame with rubber brake castors.
- Electric backrest and knee section
- FS577 Overbed lifting pole and handle.
- Single motor elevating mechanism

- Options: mattresses.
- Package size: 205cmX96cmX45cm, 99cmX36cmX12cm



» FS3230WMGF3

- ELECTRIC THREE FUNCTIONS, METAL FLAT BED PANEL
- WOODEN HEAD AND FOOT BED-BOARDS
- FOLDING ALUMINUM GUARDRAIL
- POWDER COATING STEEL FRAME
- PEDAL LOCKING CASTORS
- SIZE: 209X96X49





FS3220W

Bed with battery backup, sides and over bed pole.

- Aluminum folding sides.
- Steel epoxy coated bed frame. with rubber brake castors.
- Electric backrest and knee section.
- Fixed height bed.
- FS577 Overbed lifting pole and handle.

- Options: mattresses.
- Package Size: 205cmX96cmX45cm: 103cmX54cmX13cm

FS3220W

- ELECTRIC TWO FUNCTIONS METAL NET PANEL
- DETACHABLE ABS HEAD AND FOOT BED-BOARDS
- ALUMINUM REMOVABLE GUARDRAILS
- POWDER COATING STEEL FRAME
- PEDAL LOCKING CASTORS
- SIZE: 209X102X28

