

## Electric Obstetric Bed



B-48

### Obstetric Electric Bed

BOX SIZE 1700 X 960 X 800 mm

estimated weight 120 kg

SIZE L 2000/1450 X W 900 /660 H 620 / 920 mm

Reverse trendelenburg 0-12 degrees.

ABS Head board.

ABS Side rails.

Motors from Taiwan. Ti Motion.

Adjustable leg supports.

Hand grips.

Large bucket.

Obstetric electric bed B-48 Size:

2000/1450mm(L)\*900/660mm(W)\*620/920mm(H)

### Main feature:

1. Epoxy coated steel bed frame with ABS engineering plastic side and head board.
2. Telescopic working table with linear bearing, telescopic length : 550mm;
3. Three functions which are operated by electric motor.
4. Back rest:0 - 70±5° 2. High-low: 300±30mm.
5. Trendelenburg & reverse Trendelenburg 12 ± 2°.
6. Equiped with an 8cm thickness high density sponge mattress covered with artificial leather.
7. Adjustable leg support and sliding grip holder.
8. Four diameter-125mm castors, which are with central brakes

## Manual Obstetric Bed



B-45

### **B-45 Multi-function obstetric bed**

Box size 1400 x 890 x 910 mm

Weight + 150kg estimated

Size: 1800(L)×600(W)×650/895(H)mm

Multi-purpose obstetric table is designed for delivery,

abortion and gynecological examination.

It is characterized by easy operation,economy, safety and reliability. It is composed of bed surface, bed stead and bed base.

The bed surface includes back section, seat section and leg section.

The up-down movement of back section and the forward and backward tilt of bed surface are controlled by operating hand wheel,

which makes the doctors and patients feel comfortable.

Main specification parameters:

Length and width of table:

1800×600mm Height adjustment range: (Min)

650×(Max) 895mm

Size of sections:

Back section: 600×800mm;

Seat section: 600×400mm;

Leg section: 540×500mmMax.

degree of backrest adjustment: (up)>75°;

(down)>6°Max.

degree of forward and backward tilt of bed

surface: (forward)>10° ;(backward)>25

R

## Obstetric Examination Bed



B-43

### Epoxy coating obstetric bed

Size: 700mm(L)\*620mm(W)\*815mm(H)

Packing size:

1040mm\*900mm\*320mm (carton)

### Main features:

1. Epoxy coated frame, high quality sponge top with artificial leather cover;
2. Three sections, back rest (adjusted by gas spring), knee rest function (controlled by gear mechanic system);
3. Double leg supporting device, with one S.S. bowl;
4. Base with foot pedal, foot with anti-noise rubber

### Stationary gynaecological couch.

Epoxy powder coated tubular steel frame (32mm round tubing) mounted on 4 rubber feet.

The couch top comprises of 50mm thick, high density foam bonded to 16mm particle board and upholstered in a PVC material. Backrest, 609mm long, and leg section, also 609mm long manually operated by simple ratchets.

Fittings for lithotomy poles are provided.

Lithotomy poles included

(Dimensions; width 610mm; length 1830mm; height 762mm)



BD 881