

BABY AEROSOL NEBULISER



CODE: CNEBBABY

FEATURES:

- Compact and handy ultrasonic aerosol delivery system ideal for paediatric use
- Ideal for use on babies thank to the minimized noise level and fastness in nebulizing
- Includes an animal-shaped mask
- Light weight, easy to operate, with integrated mouthpiece and a support base for use on the table
- Can be used at "hand-held" or "desk-type"
- Automatic shut-off system after ten minutes
- Three different air flow settings to adjust treatment efficiency

SPECIFICATIONS:

| MODEL : | BABY |
|------------------------------------|--|
| Absorbed Power | 10W |
| Ultrasonic Frequency | 2.5 MHz + 200kHz (-100 kHz) |
| Size of Nebulised Particles | 1 – 10m m (measured with laser diffraction method with system SYMPATEC HELOS KF) |
| Nebulisation | Max 0.70 ml/min. – Min 0.20 ml/min (with 2ml NaCl 0.9%) |
| Functioning | 30 min. ON / 30 min. OFF |
| Auto Power Off | 11 min. (± 2 min) |
| Medication Capacity | 7.5ml |
| Maximum Sound Level | 40dB (A) |
| Weight | 0.220 Kg |
| Dimension | 60 x 117 x 164 (h) mm |