

FLO EOLO AEROSOL NEBULISER



CODE: NEBEOLO

FEATURES:

- Eolo is an electrical feed compressor with aerosol therapy nebuliser
- Designed for easy transport and handling and its recommended for atomising antibiotics and bronchodilator drugs
- The high thermal insulated plastic body complies with European Safety Standards
- Designed for non continuous use 30 min. ON / 30 min OFF

SPECIFICATIONS:

MODEL :	FLO EOLO
Power Feed	230V – 50 Hz
Power Consumption	170VA
Fuse	F1 x 1.6 A 250V
Max Pressure	2.50kPa
Operating Pressure	1.1 bar – 16p110kPa
Max Air Flow	14l/min
Operating Air Flow	5 l/min
Duty Cycle	30min ON – 30 min OFF
Noise Level	Approx. 55 Dba
Weight	1.65kg
Dimension	148mm x 124mm x 223mm
Quality Approval	CE