

# FLO MIKO AEROSOL NEBULISER



**CODE: NEBMIKO**

## **FEATURES:**

- Miko is an electrical feed compressor with aerosol therapy nebulizer
- Designed for easy transport and handling and it's recommended for atomizing antibiotics and bronchodilator drugs.
- The high thermal insulated plastic body complies with European Safety Standards
- Designed for non continuous use 20. min. ON / 40 min. OFF

## **SPECIFICATIONS:**

<b>MODEL :</b>	<b>FLO MIKO</b>
<b>Power Feed</b>	230V – 50 Hz
<b>Power Consumption</b>	170VA
<b>Fuse</b>	F1 x 1.6 A 250V
<b>Max Pressure</b>	2.2 bar – 32psi - 220kPa
<b>Operating Pressure</b>	0.95 bar – 13.8psi – 95kPa
<b>Max Air Flow</b>	13l/min
<b>Operating Air Flow</b>	4.8 l/min
<b>Duty Cycle</b>	20min ON – 40 min OFF
<b>Noise Level</b>	Approx. 56 dBA
<b>Weight</b>	1.45kg
<b>Dimension</b>	148mm x 114mm x 171mm
<b>Quality Approval</b>	CE