

OXYGEN CONCENTRATOR V3 – PURITY ALARM



<u>SPECIFICATIONS:</u>	
Oxygen flow(L/min):	1 to 3L/min
Oxygen concentration (%):	93% ±3%
Oxygen delivery pressure(Mpa):	0.05 to 0.08
Noise dB(A):	≤40Db
Size	390mm(L) x 337mm(w) x 620mm(H)
Net weight(kg):	21kg
Power:	AC 220V±10% 50Hz±1Hz
Power consumption:	≤350W
Characteristic:	Big LCD Display, accumulated working hours and the present working hours. Timing Function, Bacterial Filter, Purity alarm when purity is below 82%. Over Heat: High/Low Pressure, Low Oxygen Flow, Power Failure, Control Valve Failure
Quality Approval:	CE, ISO