

# OXYGEN CONCENTRATOR V5 SERIES 5LITRE



**CODE: OXYC5**

## **SPECIFICATIONS:**

<b>Oxygen flow(L/min):</b>	1 to 5
<b>Oxygen concentration (%):</b>	93 + 3
<b>Oxygen delivery pressure(Mpa):</b>	0.05 to 0.08
<b>Noise dB(A):</b>	≤48
<b>Atomizing pellet(m1/min):</b>	≤5μ more than 90% (Nebulizer type)
<b>Size(mm):</b>	390mm x 337mm x 620
<b>Net weight(kg):</b>	24kg
<b>Power:</b>	AC 220, 50Hz
<b>Power consumption:</b>	≤480 W
<b>Alarm Range:</b>	Overheat, high & low pressure, low flow rate, power failure, control valve failure
<b>Optional Functions:</b>	Low Purity Alarm, Nebuliser