

# OXYGEN CONCENTRATOR PORTABLE M1



**Code:**

## **Detailed Product Description**

MI portable oxygen concentrator provides continuous oxygen output.

AC or DC powered.

### **ITEM M1:**

#### **Performance:**

Power Consumption: 90W

Working Voltage: 220V±10% 50HZ±1

Car power – uses car lighter power supply (fuse max current over 20A)

Purity: 0-1L/min: 90%; 1-5L/min: 30%-90%% (Purity change according to the adjustment of flow)

Noise level: 40dB and lower

Weight: 6KG

Dimension: 324Lx199Wx320H MM

#### **Feature:**

AC/DC powered, Ideal for use in car or home.

One liter to Five liter adjust by simple switch press. Continue Oxygen Flow.

Total working hours display

Timing control function

Portable and modern design

***BACK TO OXYGEN MAIN MENU***