

VITAL SIGNS MONITOR CMS5000



OPTIONS AVAILABLE: WITH TEMPERATURE FUNCTION

CODE : PCM5000/PCM5000T

FEATURES:

- Compact and portable, allowing for uninterrupted monitoring with built-in battery
- High bright LEDS Display of NIBP, SPO2 and pulse rate
- High resolution colour LCD for trend tabular and SPO2 waveforms display
- Manual / auto / continuous measurement of NIBP
- Adjustable audible and visual alarms
- Up to 2000 groups NIBP data storage or SPO2 data up to 10 hours
- Suitable for Adult, pediatric and neonatal patient

SPECIFICATIONS:

MODEL:	CMS5000
Display	7 Segment LEDES for systolic pressure, diastolic pressure, SPO2 and pulse rate 2.4" colour LCD for waveforms, bargraphic, trend and system menu
NIBP Operation Modes: Auto Measure time: Measurement range:	Oscillometry Manual / Automatic 1 ~ 250 minutes adjustable, interval 5 minutes 10 ~ 270mmHg
Overpressure protection: Resolution: Alarm:	Yes 1mmHg Systolic, Diastolic, Mean, PR, SPO2
SP02 Measurement range: Resolution: Accuracy: Alarm range:	0 ~ 100% 1% 70 ~ 100% $\pm 2\%$, 0 ~ 69% unspecified 0 ~ 100%
Pulse Rate Range: Resolution: Accuracy:	25 ~ 250bpm 1bpm ± 2 bpm or $\pm 2\%$ (select larger)
Option:	ECG, Thermal Printer
Dimensions:	200mm x 230mm x 75mm
Weight:	Approx. 1 Kg
Power Supply:	AC100 ~ 250V 50/60Hz