

VITAL SIGNS MONITOR CMS6800



CODE :

FEATURES:

- 8" TFT, Led backlight, 800 x 480 pixels
- Visual and Audio Alarms, adjustable
- Patient Management, Name and age and ID
- Network with Central Station Software
- Audio volume adjustable form 1 ~ 8 level
- Multi-language options: Germany, France, Italian, Turkey, Etc.
- Built-in removable and rechargeable lithium battery
- UI with Both Keypad and the Touch Screen operation
- Memory: Built-in memory or miniSD Memory card. Store more than 1000 pieces of archive
- Powerful data storage and management, Trend graph, Trend table, Etc.
- Arrhythmia Analysis and S-T segment analysis
- Accuracy: $\pm 2\%$ (70% ~ 100%), 0% ~ 69% Unspecified

SPECIFICATIONS:

MODEL:	CMS6800
ECG	Lead mode: 3-lead or 5-lead HR Range: 15 ~ 300 bpm Accuracy: ± 2 bpm or $\pm 2\%$, which is greater
NIBP	Method: Oscillometry Mode: Manual/Auto/Continuous Measurement Range: 25 ~ 260mmHg
SPO2	Measurement Range: 0 ~ 100% Accuracy: $\pm 2\%$ (70% ~ 100%) 0% ~ 69% unspecified
Pulse Rate	Measurement Range: 25 ~ 250bpm Accuracy: ± 3 bpm
Temperature	Standard with Unit
Dimensions	210 x 130 x 210mm
Net Weight	2.3kg