

HANDHELD PATIENT MONITOR NIBP/SPO2 – PM50



Features

- Compact and portable, easy to use
- Suitable for adult, pediatric and neonatal patient.
- Large color LCD display of NIBP, SpO2 and Pulse rate
- Adjustable audible and visual alarms
- Real-time ambulatory monitoring. 24 hours storage and review of data.
- All data can be transferred to a PC through Windows-based software
For review and print.

SpO2

Measurement Range: 35-100%

Resolution: 1%

Accuracy: 70-100% \pm 2% 0-69% unspecified

Alarm range: 0-100%

Pulse Rate

Range: 0-254bpm

Resolution: 1bpm

Accuracy: \pm 1bpm

Alarm range: 0-254bpm

Specifications:

Display: Color LCD

Power: Two "AA", 1.5V Alkali Battery

Dimension: 120x80x32mm

Weight: <350g(Including Batteries)

NIBP

Method: Oscillometry

Mode: Manual/Automatic/continuous

Measurement Range: 25-260mmHg(Max: 280mmHg)

Automatic Measurement Interval: 5, 10, 15, 30, 45, 60, 90 minutes

Resolution: 1mmHg

Overpressure Protection: 300mmHg

Alarm: Systolic, Diastolic, mean Alarm range: 10-270mmHg