


PULSE OXIMETER HANDHELD CMS-60A



FEATURES:

- Adult Probe supplied
- Can be used in: hospital, home, community medical treatment, sports healthcare, etc.
- Can measure SPO2 and PR accurately
- SPO2 and PR display, PR waveform and bargraph Display
- Battery voltage indicator
- Friendly Operation Menu for The Function Setting
- Backlight
- Pulse sound function
- Alarm function; can set upper and lower limits.
- Interface to Computer, software can analyze, store and print the data. **(optional)**
- Function switch used to activate measurement and turn off power; function key for menu setting
- Low power consumption



CODE: POX60A

SPECIFICATIONS:

MODEL:	CMS-60A
DIMENSIONS:	109(L) X 95(W) X 40(H) MM
WEIGHT:	± 500G
DISPLAY FORMAT:	2.7" LCD display
RESOLUTION:	160*96
SPO2:	35~99%, Pulse Ratio: 30 ~ 250BPM
ACCURACY:	± 2% (70% - 99%) Unspecified (<70%) for SPO2, ±2BPM OR ±2% (Select larger) for pulse ratio
POWER:	1.5V (AA size) alkaline batteries x 2
WORKING TIME:	>30 Hours continuously
QUALITY APPROVAL	CE