

TFT COLOUR PULSE OXIMETER CMS-60C

FEATURES:

- New Look, delicate designed
- Simple operation for fast, easy test
- Superior accuracy and durability
- Shut off automatically
- Wider viewing angle with large colour LCD display
- Portable, Fashionable, Selectable colour
- Lower voltage alarm display
- Simple operation, friendly menu interface
- Flash Memory with review history data function (option)
- Interference resistance capacity against ambient light
- Can be used in Hospitals, Emergency services, Pre & Post operation, Transport, Procedure room, Nursing homes. Urgent care, sleep evaluation, Physician's office



CODE: POX60C-CS

SPECIFICATIONS:

MODEL:	CMS-60C
DISPLAY FORMAT:	LARGE TFT COLOUR LCE DISPLAY: 160 X 128
PROBES:	ADULT, NEONATE AND INFANT
POWER:	1000mAh Lithium battery, AC charger
SPO2:	35% - 99% PULSE RATIO: 30-250BPM
RESOLUTION:	±1% FOR SPO2, ±1 BPM FOR PULSE RATIO
ACCURACY:	± 2% (70%-99%), UNSPECIFIED (35%-70%), ±BPM FOR PULSE RATIO
PLETHYSMOGRAM DISPLAY	BATTERY INDICATOR
ALARM SETTING FOR SPO2 & PULSE	
QUALITY APPROVAL	CE