

PULSE OXIMETER HANDHELD CMS-60D



FEATURES:

- Routine Check Mode & Continuous Monitoring Mode
- Sync with PC based Software (WinXP and WinVista)
- Display: 1.8" Colour TFT
- Resolution: 160 x 96
- Display Mode: Display Direction adjustable
- Indicator: Power indicator, Alarm Sound, Pulse tone
- Interface: Socket for connecting SPO2 sensors and Socket for communication cables
- Power Supply: 2 AA size Batteries DC 3 V, <60MA
- Record: Data can recorded up to 24 hours
- Alarm: Adjustable High and Low limits

CODE: POX50D

SPECIFICATIONS:

MODEL:	CMS-60D
DIMENSIONS:	92MM (W) X 82MM (H) X 22mm (D)
WEIGHT:	136g
DISPLAY FORMAT:	1.8" colour TFT
RESOLUTION:	160 x 96
OPERATION ENVIRONMENT:	
TEMPERATURE:	0°C - 40 °C
HUMIDITY:	15% - 95%
STORAGE ENVIRONMENT:	
TEMPERATURE:	-20°C - 60 °C
HUMIDITY:	10% - 95%
PATIENT RANGE	Neonate, Pediatric and Adult
SP02:	
MEASUREMENT RANGE:	0 – 100%
RESOLUTION:	1%
ACCURACY:	(2%)70 – 100% (Adult & Pediatric) (3%) 70 – 100% (Neonate)
ALARM RANGE:	0% - 69% unspecified
REFRESHING RATE:	0% - 100% 1 s

PULSE RATE:	
MEASUREMENT RANGE:	25 – 250BPM
RESOLUTION:	1BPM
ACCURACY:	3PBM (NON MOTION)
ALARM RANGE:	25 – 250BPM
REFRESHING RATE:	1 s
QUALITY APPROVAL	CE