

# FINGERTIP PULSE OXIMETER FOR CHILDREN - CMS 50QA/B

**CODE :**



**50 QA – Black & White Screen Display**  
**50 QB – Colour Display**

The pulse oximeter can be used in measuring the pulse oxygen saturation and pulse rate through finger. The product is fit for family, hospital, community healthcare, physical care in sports

**FEATURES :**

- Special design for children
- Comfortable for wearing
- Low power and finger-out alarm
- Rechargeable button lithium battery
- Automatic power-off when finger is not detected

**SPECIFICATIONS :**

MODEL:	FINGERTIP PULSE OXIMETER CMS 50QA/B
Overall Dimensions	46(L) X 40(W) X 29(H) mm
Weight	± 35g (with 3.6V rechargeable lithium battery)
Power	3.6V DC~4.2VDC
Display Mode	LCD Display, SPO2, Pulse rate, Pulse rate, Pulse bar
SPO2 Measuring Range	0% - 100%
Accuracy	70 – 100%, ±,-2%, 0 ~ 69% unspecified
Resolution	SPO2: 1% Pulse Rate: 1bpm
Alarm Range	0 ~ 100% (merely applied to pulse oximeter with alarm function)
Pulse Rate	Measurement Range: 25 ~ 250BPM
Resolution	1 BPM
Accuracy	± 2 bpm or ± 2% (select larger)
Alarm Range	0 ~245bpm (merely applied to pulse oximeter with alarm function)
Quality Standard	CE