

PULSE OXIMETER WITH PRINTER - CMS65A



CODE: POX65A/PR

FEATURES:

- Simple operation, friendly menu interface
- AC/DC power supply
- Large LCD Display, low battery indicator
- Contec pulse oximeter technology
- Compatible with an array of Contec Sensors
- Portable design, convenient to carry at all times
- Impressive battery life with a quick recharge
- Built-in thermal printer, high quality

SPECIFICATIONS:

MODEL:	CMS-65A
OVERALL DIMENSIONS:	190MM X 90MM X 40MM
WEIGHT:	800G
DISPLAY FORMAT:	LARGE LCE DISPLAY; 160 X 96, LED BACK LIGHT
RECORDING SYSTEMS:	THERMAL-ARRAY; PAPER: ROLL 57MM X 25M
SPO2:	35 – 99%; PULSE RATIO: 30 – 250 BPM
ACCURACY:	SPO2 $\pm 2\%$ (70 – 99%), PULSE RATIO: ± 2 BPM
SPO2 ALARM:	0 – 100%: PULSE ALARM: 30-250 BPM
PLETHYSMOGRAM DISPLAY	
REAL TIME CLOCK	BATTERY INDICATION
OPTION:	SINGLE CHANNEL ECG MONITOR (3 LEADS), HEART RATE, HR ALARM
QUALITY APPROVAL	CE