

# PULSE OXIMETER HANDHELD PM-60A



## FEATURES:

- Routine Check Mode and Continuous Monitoring Mode
- Data Graph and Trend Table Review
- Rich Analysis Report
- Perfect Mount Solution
- Convenient operation by touch screen
- Large SD card memory for more than 2000 patient cases storage
- Sync with PC based Software (WinXP and WinVista)

**CODE: POX60A**

## SPECIFICATIONS:

<b>MODEL:</b>	<b>PM-60A</b>
<b>DIMENSIONS:</b>	92MM (W) X 82MM (H) X 22mm (D)
<b>WEIGHT:</b>	136g
<b>DISPLAY FORMAT:</b>	3.5" colour TFT Display Mode: Standard Face, Waveform face Display Direction: adjustable
<b>RESOLUTION:</b>	320 X 240
<b>INDICATOR:</b>	Power indicator light, alarm sound, Pulse tone
<b>INTERFACE:</b>	One dual-purpose socket for connecting SPO2 sensors & communication
<b>CABLES</b>	
<b>POWER SUPPLY:</b>	DC 5V, <300MA Battery Li-Polymer, 2 hours for charging, 4 hours for continuous working
<b>TREND GRAPH:</b>	Resolution from 1s, 5s, 10s. Maximum time 96 hours.
<b>TREND TABLE:</b>	Resolution form 1s, 5s, 10s. Review up to 2000 items.
<b>ALARM:</b>	Adjustable High and Low Limits. Three level audible and visual alarm
<b>OPERATION ENVIRONMENT:</b>	
<b>TEMPERATURE:</b>	0°C - 40 °C
<b>HUMIDITY:</b>	15% - 95%
<b>STORAGE ENVIRONMENT:</b>	
<b>TEMPERATURE:</b>	-20°C - 60 °C
<b>HUMIDITY:</b>	10% - 95%

<b>PATIENT RANGE</b>	Neonate, Pediatric and Adult
<b>SP02:</b>	
<b>MEASUREMENT RANGE:</b>	0 – 100%
<b>RESOLUTION:</b>	1%
<b>ACCURACY:</b>	±2% (70 – 100% Adult , Pediatric) ±3% (70 – 100% Neonate)
<b>ALARM RANGE:</b>	0% - 69% unspecified
<b>REFRESHING RATE:</b>	0% - 100% 1 s
<b>PULSE RATE:</b>	
<b>MEASUREMENT RANGE:</b>	25 – 250BPM
<b>RESOLUTION:</b>	1BPM
<b>ACCURACY:</b>	±2BPM or ±2% (select Larger)
<b>ALARM RANGE:</b>	25 – 250BPM
<b>REFRESHING RATE:</b>	1 s
<b>QUALITY APPROVAL</b>	CE