

## **BODY FAT SCALE - 150KG**



**CODE: SCABFWHC**

### **Specification**

1. Principle: Bio-electrical Impedance Analysis. (BIA)
2. Workable body weight: 3 to 150 kg / 6.6 to 330 lb
3. Increments: 100g (0.1 lb), Units: kg & lb
4. Workable Height: 100 to 210 CM (3'03" to 7'00")
5. Low battery indicator
6. Ages / Sex: 10 to 80 / Male and female. Fat measuring rang: 4% to 45 %
7. Memory space for 5 persons (9 sets each)
8. Auto step-on switch / Auto switch off function / Auto calibration function
9. Battery: "AAA" size 1.5V x 4pcs batteries (Excluded)
10. Unit Weight: 2.5kg
11. Size: 290mm x 320mm x 47mm
12. Packing: 4 pcs / CTN, Mea.: 1.49', 10.0KG (N.W.), 11.7KG (G.W.) 20' Container: 700 CTNs (2,800 pcs), 40' Container: 1,400 CTNs (5,600 pcs)

### **Features:**

- To measure body fat, it is necessary for user to setting his/her data. Different input data will result in different results.
- Display Instruction: Person number, body fat / data, age, sex, fat status.
- For next measurement, the scale detects user weight firstly and then looks for user's input data in memory by obtained weight.