

# HI-CARE DIAL BABY SCALE - MAC BA



**CODE: SCABDHC**

**Features:**

- This balance is specially made for measuring the weight of a baby
- The pan may be removed to weigh babies under 20kg standing up
- This instrument is ideal for a family, hospitals or any infant care center

**Specifications:**

<b>MODEL:</b>	<b>MAC BA</b>
<b>Capacity</b>	0 – 20kg (0 - 10kg : 50g increments / 10 - 20kg : 100g increments)
<b>Overall Dimensions</b>	545mm x 331mm x 153mm
<b>Platform Size</b>	545mm x 270mm
<b>Unit Weight</b>	3.5 kg
<b>Packing Size</b>	590mm x 330mm x 200mm
<b>Quality Standard</b>	CE