

DIGITAL BABY HANGING SCALE – MS4300



CODE:

Specifications:

MODEL:	MS4300
Capacity	25KG X 10g
Overall Dimensions	9.5 x 16cm (Scale) 45 x 31cm (Sling)
Display	0.5" LCD Display
Tare Range	Full Range
Power Supply	9V Alkaline Battery
Unit of Measure	KG/LB/OZ
Stabilisation Time	1 – 2 seconds
Operating Temperature	5°C - 35°C or 13F - 95°F
Calibration	Automatic External
Gross Weight	0.265kg
Accessories	Scale x 1, Hook x 1, Baby Sling x 1
Quality Standard	CE