

ADULT DIGITAL SCALE WITH HEIGHT ROD AND BMI – MS3400



CODE: SCA3400

Features:

- Is the ideal scale for Physicians offices, weight loss clinics , health fitness centres and other medical practices.

Specifications:

MODEL:	MS3400
Capacity	300KG X 100G increments
Overall Dimensions	270(W) x 310(L) x 70mm(H)
Platform Dimensions	310 x 310mm
Weight	6.5kg
Display	2.5cm LCD Display with 4 Digits
Power Supply	Rechargeable batteries or 15v 300mA adaptor
Key Functions	ON/OFF, kg/lb, ZERO, HOLD, BMI
Auto Off	After 120 seconds of non use
Column Height:	2M
Quality Standard	CE