

# ASKIR 20 PORTABLE SURGICAL SUCTION UNIT



**CODE : SSASK20**

## **FEATURES :**

- Easily transportable from one hospital ward to another, particularly suitable for tracheotomized patients, minor surgical applications and postoperative therapy at home
- To be used for aspirating body fluids such as mucus or catarrh
- Easily portable equipment designed for non continuous use 30 min. ON / 30 min. OFF, thanks to an electronic system for managing the electricity supply
- A series of LED's on the front control panel signals whether the machine is operating
- Made of highly heat resistant, electrically insulated plastic material in conformity with the latest European safety standards
- Supplied with a polycarbonate 1-litre autoclavable vessel complete with overflow valve
- Equipped with aspiration regulator, located on the front panel

## **SPECIFICATIONS:**

<b>MODEL:</b>	<b>ASKIR 20</b>
<b>Overall Dimensions</b>	350mm(L) x 210mm(H) x 180mm(W)
<b>Capacity</b>	<b>16l</b> per minute
<b>Unit Weight</b>	3kg
<b>Voltage</b>	230V / 50Hz
<b>Absorbed Power</b>	184 VA
<b>Fuse</b>	1 x 1.6A 250V
<b>Maximum Vacuum</b>	-0.75 BAR
<b>Standard Accessories</b>	<ul style="list-style-type: none"> <li>✓ 1l autoclavable jar (max temp. 120° C)</li> <li>✓ Cover with safety float</li> <li>✓ Silicone tubes set 6 x 10 mm</li> <li>✓ Tube connectors</li> <li>✓ Aspiration probe CH20</li> <li>✓ Antibacterial Filter</li> </ul>
<b>Quality Standard</b>	CE

***BACK TO PRODUCT INDEX***

***BACK TO SURGICAL UNITS & SPARES MAIN MENU***