

HOSPIVAC 350 SURGICAL SUCTION UNIT



CODE : SSH350

FEATURES:

- Easily transportable from one hospital ward to another
- To be used for aspirating body fluids such as mucus or catarrh and/or blood
- Easily moveable by means of wheels under the unit and a handle on top allowing virtually continuous use
- Continuous use can be enabled by activating the illuminated switch.
- Made of highly heat resistant, electrically insulated plastic material in conformity with the latest European safety standards
- Supplied with two polycarbonate 2-litre autoclavable vessels complete with overflow valve
- Equipped with aspiration regulator, located on the front panel

SPECIFICATIONS:

MODEL:	HOSPIVAC 350
Overall Dimensions	460mm(L) x 600mm(H) x 420mm(W)
Capacity	60l per minute
Functioning	Continuous
Unit Weight	13kg
Voltage	230V / 50Hz
Absorbed Power	230VA
Fuse	1 x 4 A 250V
Maximum Vacuum	-0.90 BAR
Standard Accessories	<ul style="list-style-type: none"> ✓ 2 x 2l autoclavable jar (max temp. 120°C) ✓ Cover with safety float ✓ Silicone tubes set 8 x 14mm ✓ Tube connectors ✓ Probe connector ✓ Antibacterial Filter ✓ Aspiration flow regulator ✓ Jar seat ✓ Anti-static wheels – two with brakes
Optional Accessories	<ul style="list-style-type: none"> ✓ 4l autoclavable jar (max temp. 120°C)
Quality Standard	CE

BACK TO PRODUCT INDEX
BACK TO SURGICAL UNITS & SPARES MAIN MENU