

HOSPIVAC 400 SURGICAL SUCTION UNIT



CODE : SSH400

FEATURES:

- Easily transportable from one hospital ward to another
- To be used for aspirating body fluids such as mucus or catarrh and/or blood
- Easily moveable by means of wheels under the unit and a handle on top allowing virtually continuous use
- Continuous use can be enabled by activating the illuminated switch.
- Made of highly heat resistant, electrically insulated plastic material in conformity with the latest European safety standards
- Supplied with two polycarbonate 2-litre autoclavable vessels complete with overflow valve
- Equipped with aspiration regulator, located on the front panel

SPECIFICATIONS:

| | |
|-----------------------------|--|
| MODEL: | HOSPIVAC 400 |
| Overall Dimensions | 470mm(L) x 740mm(H) x 430mm(W) |
| Capacity | 90l per minute |
| Functioning | Continuous |
| Unit Weight | 20kg |
| Voltage | 230V / 50Hz |
| Absorbed Power | 300VA |
| Fuse | 1 x 4A 250V |
| Maximum Vacuum | -0.90 BAR |
| Standard Accessories | <ul style="list-style-type: none"> ✓ 2 x 2l autoclavable jar (max temp. 120°C) ✓ Cover with safety float ✓ Silicone tubes set 8 x 14mm ✓ Tube connectors ✓ Probe connector ✓ Antibacterial Filter ✓ Aspiration flow regulator ✓ Jar seat ✓ Anti-static wheels – two with brakes |
| Optional | Foot pedal |
| Quality Standard | CE |

[BACK TO PRODUCT INDEX](#)
[BACK TO SURGICAL UNITS & SPARES MAIN MENU](#)