

DIGITAL THERMOMETER – FLEXIBLE TIP



CODE : THDFLEX

Specifications:

- Temperature range: 32°C - 42°C
- Display L°C: temperature below 32°C
- Display H°C: temperature above 42°C
- Accuracy: $\pm 0.1^{\circ}\text{C}$
- Display mode: LCD (3.5)
- Power: 0.15mW (working)
- Battery: AG3 1.5V, ordinary button battery
- Alarm: thermometer beeps after temperature taken (approx 60 seconds)
- Soft Flexible tip
- Waterproof
- Environmentally friendly - containing no mercury

Packing Details:

Single Unit – 1 complete set in Cardboard Pack
(180mm x 47mm x 15mm), 40 grams

Carton – 24 Pieces (207mm x 100mm x 190mm), 1kg

Quality Approvals:

CE ; ISO ; FDA