

Name:	Infrared Digit 5-15cm Forehead Thermometer
Model:	HTD8808
Advantage:	Non-contact Testing Forehead Temperature, Safe and Clean; Three Color Backlight (Color Alarm), Slight/High Fever Alarm; 50 Memory Storages; Alarm Temperature can be Set; 0-100°C(32-212.0°F) wide measuring range, not only for body temperature, but also for temperature of room, milk, bath water, etc with Clock function.
Measuring Range:	0.0°C~100.0°C (32.0°F~212.0°F)
Accuracy:	0.0°C-31.9°C: ±1.0°C/ (32.0°F-89.4°F: ±1.8°F) 32.0°C-35.0°C: ±0.3°C/ (89.6°F-95.0°F: ±0.5°F) 35.1°C-39.0°C: ±0.2°C/ (95.2°F-102.2°F: ±0.4°F) 39.1°C-42.2°C: ±0.3°C/ (102.4°F-108.0°F: ±0.5°F) 42.3°C-100.0°C: ±1.0°C/ (108.2°F-212.0°F: ±1.8°F)
Display:	LCD Display, Resolution is 0.1°C (0.1°F)
Alarm:	Sound alarm: sound beeps will alarm high fever; Color alarm (three-color backlight): Green backlight for normal body temperature, Yellow backlight for slight fever and Red backlight for High fever.
Units:	°C/°F can be Changed.
Memory:	50 Memory Storages
Test Mode:	5-15cm Non-contact Testing Forehead Temperature
Measure Time:	1 Second
Auto Shut-off:	5 Seconds without Operation
Battery:	2 pcs AAA batteries
Battery Life	Approximate 3,000 Times, Inspect the Low Voltage of Battery at any Time.
Operating Temperatures:	5.0°C-40.0°C(41°F-104°F), Relative Humidity≤80%
Storage Temperatures:	-20.0°C~70.0°C(-4°F-158°F), Relative Humidity≤95%
Dimensions:	150 x 95 x 44mm
Weight:	115g (without Battery)
Standards:	CE,ROHS,FCC

