

# ULTRA VIOLET STERILISER 8 WATT OR 15 WATT



**CODE : UVS8/UVS15**


## FEATURES :

- Violet is a germicidal system which uses an ultraviolet germicidal lamp
- This Violet maintains sterility of instruments and small disposables after sterilization in an autoclave or similar type of sterilization unit
- Stainless steel body
- Inexpensive, reliable, safe, easy operation and maintenance

## SPECIFICATIONS :

MODEL :	Ultra Violet
Overall Dimensions	8W – 195mm x 165mm x 375mm 15W – 270mm x 215mm x 495mm
Types	DC510000 – VIOLET 8W / DC510000 – VIOLET 15W
Absorbed Power	8W / 15W
Fuse	2 x 1Amp
Optional Accessories	✓ 8 Watt Lamp / 15 Watt Lamp ✓ Grill for Germicidal Lamp
Quality Standard	CE

**8W**



**15W**

