



# PRODUCT GUIDE

Your local distributor:

**WelchAllyn®**

Advancing Frontline Care™









## WELCH ALLYN – SPECIALISING IN:

### Table of Contents

#### Diagnostic Products Guide

■	Diagnostic Sets and Instruments	4
■	Diagnostic Sets and Power Handles	5
■	Wall and Desk Systems	6
■	Diagnostic Sets	7
■	Stethoscopes and Laryngoscopes	8
■	Blood Pressure	9
■	Blood Pressure and FlexiPort Universal Cuffs	10
■	Screening	11
■	Thermometry	12
■	Headlights	13
■	Womens Health	14
■	Lighting	15
■	Vital Signs Spot	16
■	Continuous Monitoring	17
■	Cardiopulmonary ECG	18
■	Cardiopulmonary	19





## FRONTLINE CARE FOR A CENTURY.

[CLICK HERE TO VIEW VIDEO](#)

Get everything you need for frontline care with Welch Allyn.

You can't always predict how many patients you'll treat or the type of care they'll need.

Welch Allyn has the equipment you need to meet challenges such as these and give you more time to focus not only on the initial diagnosis but also ongoing care. For 100 years, we've been helping frontline care givers like you treat more patients, more efficiently. Our family of products—from handheld diagnostic tools to sophisticated patient monitors—are renowned for their simplicity and durability.

This product guide provides you with a snapshot of some of the equipment available from Welch Allyn. Discover for yourself why Welch Allyn is a name you can count on in any clinical setting.

- Focused exclusively on frontline care
- Worldwide sales, service and support
- Simple, intuitive devices that take very little of your time to learn
- Renowned for durability and reliability
- Trusted by thousands of healthcare professionals every day







23820 / 23810

## MACROVIEW™ OTOSCOPE

- Delivers almost twice the field of view at higher magnification, making it easier to view
- The ability to adjust focus for variable ear canal length or far sighted eyes
- Fibre optics produce cool light with no reflections or obstructions

[CLICK HERE TO VIEW VIDEO](#)



20000

## 3.5V DIAGNOSTIC OTOSCOPE

- Halogen HPX™ lamp provides 30% more light output than standard halogen lamps for true tissue colour
- Fibre optics project cool light to the distal tip with no reflections or obstructions for superior views of the tympanic membrane
- Sealed system allows for pneumatic otoscopy and accurate diagnosis of otitis media

[CLICK HERE TO VIEW VIDEO](#)



23920

## DIGITAL MACROVIEW™ OTOSCOPE

- Improves communication which can lead to better patient compliance
- Allows storage of digital images in patient records for enhanced documentation of exam findings
- Provides the ability to share images with other providers as visual descriptive reference
- Gives an excellent training aid for educational purposes
- Offers the ability to view small details of tubes and hearing aids

[CLICK HERE TO VIEW VIDEO](#)



## PANOPTIC® OPHTHALMOSCOPE

- 5 X larger view of the fundus vs standard ophthalmoscopes in an undilated eye
- 25° field of view vs. the standard 5° field of view
- Greater working distance improves comfort for both practitioner and patient
- Patented Glare Extinguishment System prevents interference from unwanted glare and reflections
- Take photos with your iPhone



CLICK HERE TO VIEW VIDEO

CLICK HERE TO VIEW VIDEO



## 3.5 V COAXIAL OPHTHALMOSCOPE

- Coaxial optics produce a shadow-free spot, easier entry into undilated pupils, and a larger field of view versus standard ophthalmoscopes.
- Detect corneal abrasions with cobalt blue filter
- Polarizing filter virtually eliminates corneal reflection
- 18 unique aperture/filter combinations for greater versatility
- Red free filter may be used with any aperture
- 28 focusing lenses with a range of -25 to +40 diopters

CLICK HERE TO VIEW VIDEO



## LED REPLACEMENT BULBS

- Long life backed by a 7-year warranty (date stamped)
- SureColor™ LED technology enabling reduction of intensity to alleviate glare or enhance patient comfort while preserving the colour and clarity of the view
- High energy efficiency resulting in 3X the examination time with each charge of a Welch Allyn rechargeable handle

CLICK HERE TO VIEW VIDEO



- 03100 - LED
- 03800 - LED
- 04900 - LED
- 06500 - LED

## 3.5 V NICKEL-CADMIUM RECHARGEABLE/ CONVERTIBLE HANDLE

- Smooth, non-abrasive knurled finish ensures a no-slip grip
- Chrome-plated brass construction for durability and reliable long life
- Re-chargeable or can be converted for use with standard C-cell battery handles (model 71020-C)



## 3.5 V LITHIUM ION RECHARGEABLE HANDLE

- Half the weight of the traditional 3.5V handles
- Fits all Welch Allyn 3.5V instrument heads
- 120 minutes of on-time as opposed to the 60 minutes for standard handles
- Lithium ion battery technology that has no memory effect





## GS™ 777 WALL TRANSFORMER

- Help reduce facility cost with a more energy efficient product than the 767 Wall Transformer
- Help improve workflow and save energy with handles that automatically turn on/off
- No on/off switch reduces contact, to help reduce the risk of cross-contamination
- Upgrade to the Welch Allyn LED lamp for greater energy efficiency and lower cost of ownership



77747



77797-M2X

## 777 INTEGRATED WALL SYSTEM

In addition to the 777 Wall System the fully integrated system also includes:

- Wall Aneroid (7670-01)
- Disposable Specula Dispenser (52401F)
- OPTIONAL EXTRA—SureTemp® Plus Electronic Thermometer
- Can be configured to each facility's needs



## GREEN SERIES™ DIAGNOSTIC DESK SYSTEM

- Desk Aneroid w/Adult Blood Pressure Cuff
- 3.5V Lithium Ion or 3.5V Ni-Cad Handle
- MacroView™ Otoscope Head w/ Throat Illuminator, Braun ThermoScan PRO 4000 Thermometer
- Space for optional second handle with Ophthalmoscope or PanOptic
- Storage and charging station in one convenient location
- Anti-theft function (Pro4000)



## PROFESSIONAL PLUS DIAGNOSTIC SET

- Coaxial Ophthalmoscope (11720)
- Diagnostic Oscope (25020)
- C-cell battery handle (71020-B)
- Hard case (05259)

[CLICK HERE TO VIEW VIDEO](#)



97250-BI

## HOSPITAL CLASSIC DIAGNOSTIC SET

- 3.5 V Coaxial Ophthalmoscope (11720)
- 3.5 V Oscope—HPX™ with Illuminator (20000)
- NiCad rechargeable/convertible handle (71062-C)
- Hard case (05259)



97202-C

## HOSPITAL ELITE DIAGNOSTIC SET

- 3.5 V Coaxial Ophthalmoscope (11720)
- 3.5 V Diagnostic Oscope with HPX™ illumination and illuminator (20000)
- Lithium ion rechargeable Smart handle (71902)
- Hard case (05259)



97202-VSM

## NEW POCKET LED SET

- Welch Allyn SureColor™ LED lamps
- Most durable Pocket set on the market today
- Lightest-weight Welch Allyn otoscope and ophthalmoscope
- Fiber optics provide cool light with no reflections or obstructions
- All the apertures and filters as in our full-size ophthalmoscopes



92880 - BLK or WHT or BLU or PUR



92871 - BLK or WHT or BLU or PUR

[CLICK HERE TO VIEW VIDEO](#)





5079-125S

## HARVEY™ ELITE® STETHOSCOPE

- Double head (bell and flat diaphragm) high quality stainless steel chest piece with non-chill rim
- Rotatable stainless steel binaurals and interchangeable comfort sealing ear tips
- Dual channel PVC tubing provides separate sound pathways to each ear for unrivalled accuracy
- Triple leaf binaural spring encased in polyurethane allows size adjustment while preventing breaking
- Paediatric chest pieces available
- Latex free



5079-135

## PROFESSIONAL STETHOSCOPE

- Meets all diagnostic needs, from detecting Korotkoff sounds of blood pressure to analyzing subtle heart murmurs
- Comfort sealing eartips depress and conform to ear canal for soft yet airtight fit
- Latex free
- 5 year warranty

## FIBRE OPTIC LARYNGOSCOPE HANDLES

- Patented HPX™ technology blends halogen and xenon gases under pressure to produce a bright white light for true tissue colour
- Up to 40% lighter for increased manoeuvrability
- C-Cell Handle, AA Handle, Stubby Handle
- **NEW** LED bulbs



60813

60814

60815

## FIBRE OPTIC LARYNGOSCOPE BLADES

- Single piece stainless steel construction that ensures blade integrity
- Damaged light pathways can be repaired at a fraction of the cost of replacement blades
- Miller blades—size 00 to 4,
- English Macintosh—size 1 to 4
- Macintosh blades—size 1 to 4



Macintosh



Emac



Miller



## DURASHOCK™ HANDHELD OR TRIGGER MANOMETERS

- Welch Allyn Durashock technology—Shock resistant for accuracy and lower cost over time—gear free—remain in calibration for longer
- Will survive a 76cm drop and remain in calibration—AAMI shock resistant specification
- Jewel movement for long life
- Laser engraved dial for accuracy
- 5 (Bronze), 10 (Silver) and 15 (Gold) year warranty

[CLICK HERE TO VIEW VIDEO](#)



DS66 GOLD RANGE



DS-5401-189 BRONZE



DS-5501-189 SILVER

## 767 SERIES - WALL, MOBILE AND DESK TOP

- Certified accuracy to  $\pm 3$ mmHg
- Mercury and latex free
- Ruby jewelled movement for long life
- Large, easy to read dial face—15cm
- Built-in cuff storage basket
- Wall mount option
- Mobile stand option—adjustable height and 360° swivel
- Desk stand option
- Lifetime calibration warranty

7670-01



## NEW PROBP 2400

- Angled desktop design with an easy user interface and backlit LCD display.
- Auto inflate mode senses proper maximum inflation pressure for improved patient comfort.
- **MODES:**
  - Single Reading Mode for routine blood pressure screening.
  - Three-Reading Average Mode helps eliminate “white-coat” effect.
  - Manual Mode allows traditional auscultatory blood pressure method. ProBP 2400 acts as digital aneroid gauge. Auto inflates then deflates at 3 mmHg/min.
- 4 Cuff sizes from age 12 and above.
- Folding Handle



2400

### MOUNTING OPTIONS:



140D-WMB



2400-MS



[CLICK HERE TO VIEW VIDEO](#)

## CONNEX® PROBP™ 3400 DIGITAL BLOOD PRESSURE DEVICE

- Small handheld design
- Lithium ion battery delivers 100 readings on a single charge
- Compatible with blood pressure cuff sizes ranging from small child (12cm) to thigh (55cm)

### MOUNTING OPTIONS:



34XFHT-7

[CLICK HERE TO VIEW VIDEO](#)

## FLEXIPORT® BLOOD PRESSURE CUFFS

- Easy one handed attachment/detachment of tubing from FlexiPort makes changing of cuffs easier and quicker
- Standardizing can reduce your cuff part numbers by up to 60%
- Rotatable port reduces stress to cuff and tubing, increasing durability

### The right cuff. The right tubes. At the right time.

The unique FlexiPort fitting enables true cuff standardization throughout your facility for the first time. With FlexiPort, finding the right cuff is made simpler than ever before by eliminating annoying hassles and frustrations of tube and connector configurations.

[CLICK HERE TO VIEW VIDEO](#)

## FLEXIPORT® ECOCUFF™

- Improve revenue and Profitability with cuff Standardisation
- Single-patient-use FlexiPort® EcoCuff™
- FlexiPort® single point connection technology for controlling cross-contamination
- Cover a full range of sizes and meet applicable AAM standards





[CLICK HERE TO VIEW VIDEO](#)

## NEW SPOT™ VISION SCREENER

- Spot supports AAP vision screening guidelines for early detection of amblyoptic risk factors
- Easy to use
- Lights and sounds designed to help to engage children during the screening
- Spot's high capture rate produces meaningful results
- WiFi-enabled for easy wireless printing and exporting of data



VS100-7

## NEW RETeval-DR™

- Screening for early detection of Diabetic Retinopathy and clinically significant macular edema (CSME) screening
- Hand-held and battery powered
- Instant results displayed on screen
- Built-in pupillometry adjusts flash and background intensity, enabling testing without pupil dilation
- Single electrode array is placed on skin for each eye (Patent Pending)
- IR camera allows visualization of eye during test
- EMR interface available
- Reports viewable on the device or on a PC using the USB interface in the base unit



RETeval-001

## NEW OAE HEARING SCREENER

- Welch Allyn OAE Hearing Screener offers objective hearing screening technology, and tests at key speech frequencies (2, 3, 4, 5 kHz)- ideal when testing newborns, infants, toddlers, preschool and school-age children
- Low-Cost Disposable Ear Tips from 3 mm to 15 mm.
- Durable Probe: Designed with HDMI connector to help eliminate connector pin breakage.
- Wireless Printer: Bluetooth thermal printer provides graphic and tabular reports in seconds with a single button press.



39400





01690-410

05031-101

## SURETEMP® PLUS THERMOMETER

- Interchangeable, removable probe wells (rectal and oral) reduce the risk of cross contamination
- Includes a 60 second pulse timer for counting pulse and respiration
- Mechanical or electronic security options to help prevent theft
- Takes approximately 5000 readings on three AA batteries

[CLICK HERE TO VIEW VIDEO](#)



04000-600

04000-800

## BRAUN® THERMOSCAN® PRO 4000 EAR THERMOMETER

- Improved accuracy and repeatability through a preheated probe tip and sensor
- Probe cover detection system ensures precise infrared heat transfer for accurate temperatures
- Automatic probe cover eject button for quick and easy disposal of used



06000-200

## NEW BRAUN THERMOSCAN® PRO 6000 EAR THERMOMETER

- Advanced measurement technology, allows quick and accurate temperature readings
- Pre-warmed sensor helps support accurate measurements
- 60-second pulse timer, memory recall button, C/F scale conversion button
- Pre-warmed probe tip enhances patient comfort during reading
- Probe covers are soft-walled and disposable

[CLICK HERE TO VIEW VIDEO](#)



105801

## NEW CARE TEMP™ TOUCH FREE THERMOMETER

- Infrared sensor reads temperature 4-6 cm from the center of forehead
- No patient skin contact required, helps minimize the risk of cross-contamination
- 2 year device warranty
- 60 second pulse timer, memory recall button, C/F scale conversion button
- Measurement feedback system detects proper operating distance for accuracy
- Large, backlit LCD display that's intuitive and easy to read
- Takes approximately 3,000 readings on two AAA batteries



## LUMIVIEW™ PORTABLE BINOCULAR MICROSCOPE

- Unique 3-D optics or viewing into small cavities
- Easy adjustable headband
- Cool halogen light
- Portable power source with charger

[CLICK HERE TO VIEW VIDEO](#)



20502 H

## SOLID STATE PROCEDURE HEADLIGHT

- Lightweight
- Colour temperature do 5500 degree Kelvin
- Estimate life of up to 10,000 hours continuous burn
- Direct power or portable power source



49022

[CLICK HERE TO VIEW VIDEO](#)

## SOLID STATE PORTABLE HEADLIGHT

- Lightweight
- Colour temperature do 5500 degree Kelvin
- Estimate life of up to 10,000 hours continuous burn
- Direct power or portable power source
- Battery run time: Approximately 90 minutes of continuous run-time on a fully charged battery
- Solid State Portable Headlight with Direct Power Source / Charger
- Available in Soft Terrycloth or Rigid Headbands



46072



46072R







79912

## KLEENSPEC ILLUMINATION SYSTEM

- Advanced LED light sources provides improved visualization of the exam area
- No cord to get in the way during procedures
- Automatically illuminates upon insertion of the illuminator and once its removed the light goes off

[CLICK HERE TO VIEW VIDEO](#)



## KLEENSPEC VAGINAL SPECULA

- Newly designed light pipe provides a homogenous spot and eliminates glare back
- Smooth moulded acrylic will not pinch or bind tissue for greater patient comfort during the exam
- Available in sizes extra small, small, medium and large





[CLICK HERE TO VIEW VIDEO](#)

## LIGHTING - GREEN SERIES RANGE

- 5,500° K white light-emitting diodes (LEDs) for brighter, whiter light
- Truer tissue rendition and uniform spot (no dark or hot spots)
- Lower cost of ownership—no bulb replacement and low energy consumption
- Touchless on/off control



LED VS. HALOGEN

### GS EXAM LIGHT IV

- Perfect for OB/GYN office or for specialities that require focused light
- 3X the typical light output – highly flexible light pipe
- Table, wall or mobile stand
- 5 year limited warranty



48817

### GS 300 GENERAL EXAM LIGHT

- Ideal for primary care facilities, emergency rooms, outpatient facilities and dermatology
- Table, wall or mobile stand
- 5 year limited warranty



44457

### GS 600 MINOR PROCEDURE LIGHT

- Perfect for minor procedures and exams in any healthcare environment
- Features 3 LED's and a wider head with single handle
- Touchless dimming control
- Table, wall or mobile stand
- 5 year limited warranty



44607

### GS 900 PROCEDURE LIGHT

- Ideal for procedures and exams in any hospital or healthcare environment
- Features 6 LED's and a wider head with double handles
- Touchless dimming control
- Table, wall, mobile stand or ceiling mount
- 5 year limited warranty



44900-C





71WT-7  
74CT-7  
75WT-7

## NEW CONNEX SPOT MONITOR

[CLICK HERE TO VIEW VIDEO](#)

- Help improve accuracy and performance for neonate through adult patients
- Blood pressure averaging to help ensure accurate hypertension diagnosis
- Choice of leading SpO2 technologies: Masimo®, Nellcor® or Nonin®
- Mobile stand option incl. work surface w/integrated battery for extended usage
- SureBP technology can provide BP in only 15-seconds - 3X faster than competitors
- Connex Scoring App on-device calculator for up to 3 custom protocols like NEWS, MEWS
- Partnerships with more than 90 leading EMRs including Epic and Cerner



45NTO-E2

## SPOT VITAL SIGNS® LXI

- SureBP™ technology measures Blood pressure in 15 seconds
- Choice of Suretemp Plus or Braun thermometry
- 50 reading memory allows for retrieval of previous data
- Configurations to include any of the following: NIBP, Temperature and SPO2

## SPOT VITAL SIGNS®

- Measures NIBP, pulse rate, temperature and SPO2 in 30 seconds
- Lightweight and portable with a variety of mounting options
- All configurations include NIBP, pulse rate and MAP
- Optional Masimo® or Nellcor® SPO2
- Optional SureTemp® thermometry



420NTB-E2

[CLICK HERE TO VIEW VIDEO](#)





## PROPAQ® LT

- Motion Tolerant Monitoring
- Engineered to withstand shock up to 75g and a drop of 2 meters
- Suitable for both bedside, transport and ambulatory
- Long battery life—between 8 and 24 hours
- Wireless connectivity option
- Large colour display option
- 3-5 lead ECG
- Respiration
- SPO2
- NIBP



## CONNEX® INTEGRATED WALL SYSTEM

- Combines advanced spot check and monitoring capabilities with an otoscope, ophthalmoscope and thermometer all in one convenient location
- For connected facilities, vitals can be electronically sent directly to the EMR/EHR system
- SureBP-fast, 15-second NIBP and pulse rate
- Braun ThermoScan® PRO 4000 ear thermometry
- Masimo® or Nellcor® pulse oximetry
- Adult, pediatric, and neonatal blood pressure modes
- NIBP interval monitoring mode, including customizable interval programs
- Programmable visual and audible alarms for BP, SpO2, temperature and pulse
- Integrated 3.5 V handles for physical assessment tools
- Optional Barcode Scanning for Patient Identification
- Optional wireless integration to Seca scales for digital height and weight measurement



## CONNEX® VITAL SIGNS MONITOR

- Specially designed for non-critical care settings
- Acts as 3 devices in one: spot check, monitoring and triage
- Vivid 8.9 inch colour touch screen—facilitates cleaning and prevention of infections
- Smart Lithium Ion battery for longer "on" time
- Configurations to include any of the following: NIBP, Temperature, SPO2 and Printer
- Accessory and cable management stand with unique cable routing system all components have their own compartment for quick ID and access
- Optional Early Warning Scoring System
- Optional Barcode Scanning for Patient Identification
- Optional Respiratory Monitoring
- Optional wireless integration to Seca scales for digital height and weight measurement



## CONNEX® CENTRAL STATION

- Unique numeric-based user interface replicates the Connex Vital Signs Monitor
- Renew vital sign and continuous monitoring in tabular or graphical format
- Visual and audible indications from central station and hallway displays
- EMR integration helps to automatically document accurate patient vitals
- Supports up to 48 patients; multiple central stations can be networked together
- Intermittent vitals and continuous monitoring combine in a single patient record
- Intuitive review of retrospective patient data and trends





[CLICK HERE TO VIEW VIDEO](#)



CP50-2EN9-BASIC  
CP50-2EN9-PLUS

### CP 50™ RESTING ECG

- Non-interpretative 12 lead ECG—interpretation software optional extra
- 110mm colour touch screen
- 100mm high resolution thermal printer
- “Lead off” and quality messages
- Lithium Ion Smart battery
- Full charge in less than 3 hours
- External printer support option
- Disposable or re-usable electrodes

[CLICK HERE TO VIEW VIDEO](#)




CP150-2EN7  
CP150A-2EN7  
INTERPRETATION

### CP 150™ ELECTROCARDIOGRAPH

[CLICK HERE TO VIEW VIDEO](#)

- 7” Colour, touch-screen display for fast and easy entry
- Instant-on feature powers the device up quickly and lets you take an ECG at a moment’s notice
- Advanced filters for optimal ECG trace quality
- Optional adult and paediatric MEANS ECG interpretation software
- Battery operation lets you easily transport your device to the patient
- Store up to 100 test results on the device or transfer to a USB memory stick
- Optional Spirometry



STRESS:  
CPS-UN-EP

RESTING  
CPR-UN-EB

### PC-BASED RESTING ELECTROCARDIOGRAPH

- Use with a single computer or enterprise wide network, with or without an EMR
- Improve exam efficiency with a 4x3 preview on screen prior to printing
- Automatic ECG measurements speed physician review and interpretation
- Help improve patient care by overlaying previous ECGs to spot differences
- Interpret and edit report on screen for patients down to 1 day old
- Easily monitor patients with intermittent arrhythmias
- Programmable report formats

[CLICK HERE TO VIEW VIDEO](#)

### PC-BASED EXERCISE STRESS ECG

- Select any of the commonly used protocols or customize your own
- Computer based—low cost stress test system
- 12 Lead ECG mode for baseline ECG’s
- Records and monitors all 12 leads simultaneously
- Reports printed on demand or at pre-defined intervals
- Arrhythmia detection is available throughout the test
- All data is stored on the hard drive
- Optional Spiro, Hoter, ABPM, Interpretation



721070  
Treadmill with Controller

BI-WA  
Ergometer



## SPIROPERFECT™ SPIROMETER

- Fast, reliable Pre-Post bronchodilator FVC, SVC and MVV testing with automatic interpretation and reversibility analysis
- View real-time Flow/Volume Loop and Volume/Time graphs
- Instantly verify quality and variability of proper test performance
- Motivate children with incentive graphic - for pediatric patient coaching
- Reduce risk of cross contamination with single-use disposable transducers
- Select single and multiple flow calibration protocols

[CLICK HERE TO VIEW VIDEO](#)



## NEW ABPM 7100 AMBULATORY BLOOD PRESSURE MONITOR

- Meets or exceeds the latest AAMI, EHS Standards and has a BHS A/A Rating
- Unlimited protocols for programming measurement periods & inflation frequencies
- Lightweight, compact and quiet for improved patient compliance
- Uses oscillometric technology
- Upgradable to add Central BP and Pulse Wave Analysis
- Auto Feedback Logic during inflation for excellent patient comfort & acceptance
- Latex-free cuffs help reduce the risk of patient allergic reactions
- Includes Welch Allyn CardioPerfect® WorkStation software



## HR-100 HOLTER RECORDER

Holter recording in an incredibly small package

- 2- or 3-channel recording
- 24-48 hours of recording time
- Pacer detection
- Wireless transfer of patient's ECG in real time to verify electrode quality and placement
- Optional QT Dispersion Software



## HR-300 HOLTER RECORDER

Powerful features for the demanding user

In addition to the features offered in the HR-100, the HR-300 offers the following:

- Display for real-time ECG preview or test programming
- Remote initialization helps to enhance workflow
- Programmable sampling rates of 200, 500 and 1,000 sps





## CONSUMABLES AND ACCESSORIES

Product Name	Part No.	Ref. Page	Product Affiliation
Thermometry	#05031-101	10	Disposable Thermometer Probe Covers for SureTemp® 678, 679, 690, 692 and 986 Models; Spot Vital Signs®LXi; Spot Vital Signs® (1,000 each)
	#04000-800	10	Disposable Thermometer Probe Covers for Braun ThermoScan® PRO 4000 (800 each)
Paper	#406021	16	CP50 Paper
	#52600	11	TM 286™ AutoTyp™ Printer Paper (1 pack of 5 rolls)
	#56100	11	MicroTyp™3 Printer Paper (5 rolls)
	#105353	16	Paper for CP 100™/ CP 150™/ CP 200™
Electrodes	#45008-0000	n/a	CP 300 Resting Tab Electrodes (box of 1,000)
	#715006	16	Multi-Functional Electrode Adapters (10 clips)
Batteries	#71032	n/a	Charger 220-240v Input Compliant with IEC 601-1
	#72300	5	Replacement NiCad Rechargeable Battery (orange) for 71000-A/71000-C
	#72200	5	Replacement NiCad Rechargeable Battery (black) for Desk or Wall Charger Handle
	#71960	5	Replacement Battery for Lithium Ion Handle
Lamps	#03000-U	n/a	Ophthalmoscope #11710
	#04900-U	4	Ophthalmoscopes #11720, 11730, 11735
	#04900-LED	4	Ophthalmoscopes #11720, 11730, 11735
	#03800-U	5	PanOptic™ Ophthalmoscopes #11810, 11820
	#03800-LED	5	PanOptic™ Ophthalmoscopes #11810, 11820
	#06500-U	4	MacroView™ Otoscope #23810, 23820
	#06500-LED	4	MacroView™ Otoscope #23810, 23820, 23920
	#03100-U	n/a	Otosopes #20000, 20200, #21700, 20200, 3.5 V Handle Adapter #73500
	#03100-LED	n/a	Otosopes #20000, 20200, #21700, 20200, 3.5 V Handle Adapter #73500
	#06000-U	4	Halogen HPX replacement lamp
	#04700-U	19	MacIntosh standard lamp blades, length 89mm or 108mm
	#04800-U	19	MacIntosh standard lamp blades, length 126mm or 160mm
	#08500	11	LumiView bulb
	#08800-U	12	4.6 V Halogen Lamp for Vaginal Specula Illumination System #78810
	#09800-U	12	Solarc® Lamp for Video Colposcope
Otoscope Specula	#24400	n/a	Renewable Spec Set (4)
	#52432-UB	5	24 Series, MacroView, Disposable KleenSpec,® 2.75 mm (850/bag, 10 bags/case)
	#52434-UB	5	24 Series, MacroView, Disposable KleenSpec,® 4.25 mm (850/bag, 10 bags/case)
	#52700	4	Slotted Tip for #23810 and #23820 MacroView Otoscope (48/box, 4 boxes/case)
	#24320	4	MacroView,™ Disposable SofSeal™ Sealing Accessory for Otoscope, Small, for 2 to 4.25 mm (80/box, 6 boxes/case)
	#24330	4	MacroView,™ Disposable SofSeal™ Sealing Accessory for Otoscope, Medium, for 3 mm, 4 mm, 5 mm (80/box, 6 boxes/case)

Welch Allyn South Africa (Pty) Ltd  
 Building 9-Pebble beach  
 Fourways Golf Park, Roos Street  
 Fourways, 2191



**WelchAllyn®**

Advancing Frontline Care™

Tel: +27 (0) 100017788 | Fax: +27 (0) 11 465 4576  
 www.welchallyn.co.za